ENGLISH

SpineVision®

Purpose:

Flex+2 implants are designed for the correction and spinal stabilization of the spinal column only in combination with pedicle screw systems developed and marketed by SpineVision: titanium PL.U.S.™, LUMIS® or U.L.I.S.™. The rods have been designed to allow the surgeon to adapt the stabilization type to the degree of degeneration.

Description:

The Flex+2 range comprises different types of rods: dynamic stabilization rods (F1-RH609T) and hybrid stabilization rods (F1-RH609T) to associate an arthrodesis with a dynamic stabilization thanks to the flexible part of the rod on one level (F1-RH1609T) or two contiguous levels (F1-RH2609T).

The metallic components are made of titanium alloy (Ti-6Al-4V) compliant with standard ASTM F136 (ISO 5832-3), and the parts used for dynamic stabilization are overmolded with biocompatible Urethane Polycarbonate (USP class VI).

WARNING: The Flex+2 rods are not to be used with screws or instruments from another manufacturer, or with pedicle systems not described in the instructions.

According to the clinical evaluation report, it can be concluded that the benefit-risk ratio is positive.

- Indications:**
- When properly used on a patient whose skeleton is mature, any flexible part of the Flex+2 rod is intended to stabilize vertebral segments over one or two adjacent stages from T10 to L5 in the following spinal pathologies:
- Moderate lumbar discopathy at Pfirmann stages II or IV
 - Early degenerative spondylolisthesis (grade I) and hyper-mobile segment without extensive facet joint osteoarthrosis
 - Lumbar stenosis by bulging of discs and of flwal ligaments, associated with a recalcitrant,
 - Massive herniated disc in a young patient with anticipated discophy
- Any rigid part of the Flex+2 rod is intended to fuse several vertebral segments from T11 to S1, directly adjacent to the dynamically stabilized levels in the following spinal pathologies:
- Degenerative discopathies
 - Facet degeneration
 - Degenerative spondylolisthesis
 - Lumbar stenosis
 - Pseudarthrosis following a previous instrumented surgery.

Contraindications:

Contraindications include:

- All cases not described in the indications.
- Treatment of thoracic spine above T10 and of cervical spine.
- Unilateral application of the Flex+2.
- Use of two rods with different references on the same level
- Use of PL.U.S.® titanium screws, LUMIS® or U.L.I.S.™ screws with a diameter below 5.5 mm at the flexible level,
- Any previous surgery at the level(s) or previous instrumented lumbar surgery at the rigid level(s) except if pseudarthrosis can be proved,
- Pathological obesity (Body mass index, BMI > 35),
- Pregnant woman should not undergo a surgery.
- Concomitant performance of facetectomy and/or laminectomy at the level(s) instrumented with the flexible part(s)
- Known or suspected allergy, metal and/or urethane polycarbonate intolerance.
- Localized infection of the surgical site, signs of inflammation, fever, systemic or local acute infection.
- All patients with insufficient tissue coverage of the surgical site.
- Radically developing arthropathy, acceleration of bone resorption, osteopenia / osteoporosis or osteomalacia highlighted by X-rays or fracture antecedents. The presence of osteoporosis is a relative contraindication, as it can limit the expected correction gain and the stability of the mechanical fastening.
- Smokers whose bone density is below a T-score of -1.5.
- Arthrosis at the flexible level(s)
- Degenerative scoliosis > 10° at affected level(s).
- Spondylolisthesis higher than grade I and spondylolisthesis by spondylolysis at the flexible level(s).
- Any other medical or surgical affection liable to compromise the success of the instrumented operation, such as, extensive deformation, tumor, vertebral fractures or serious congenital anomalies, acceleration of the erythrocyte sedimentation speed not accounted for by another disease, high leukocyte count or tendency to leukopenia.
- All cases requiring the association of different metals.
- Patient's inability to follow practitioner recommendations (Mental illness).
- Any patient who does not accept to follow the postoperative instructions.

The contraindications for these devices are similar to those for other instrumentations by vertebral spinal rods. This spinal instrumentation device is not designed or sold for applications other than those mentioned.

WARNING:

Use of the Flex+2 on the flexible part(s) of the Flex+2 implant is prohibited to avoid damaging the overmolded part. For hybrid rods, the flexible part must always be implanted at the cranial end of the assembly.

This device can only be implanted by a surgeon who is thoroughly familiar with it and with its applications, and with the instruments and the surgical technique required.

The surgeon must have a sound knowledge, not just of the medical and surgical aspects of the implant, but also of the material limits of implants in urethane polycarbonate.

The patient must be informed of the limits of the implant, particularly in terms of loads and other bodily stresses exerted on the device prior to forming of a solid arthrodesis with fused stages.

The patient must also be aware that failure to follow the postoperative instructions could lead to failure of the device, possibly requiring a new operation to remove it.

Patient selection:

- Particular attention must be paid to the patient selection criteria, to the proper setting of the implant, and to operative care, in order to minimize the stresses to which the implant is subjected. In particular, it is important to take into consideration the patient's sagittal equilibrium so as to take his pelvic inclination parameters into account and avoid excessive shear stresses on the flexible part(s).
- Only the patients meeting the criteria described in the indications can be selected.
- Patients who meet the criteria described in the contraindications must not be selected.

Potential adverse effects:

The following adverse events are inherent to any surgical operation, or associated with instrumented intervention on the spine, including in particular:

- Coming apart, warping, slipping and/or rupture of any or all of the components or instruments during surgery.
- Reaction to the implant as a foreign body, e.g. possible tumor formation, autoimmune disease, metallosis and/or healing complications.
- Decrease in curvature or correction of the spinal column or in the height of the discs.
- Infection.
- Fracture of the vertebral body due to load transfer above and below the instrumented area.
- Absence of bone setting (pseudarthrosis) at the rigid part.
- Pressure on the skin or other soft tissues from component parts possibly causing skin penetration, irritation, fibrosis, necrosis, and/or pain.
- Device breakage, possibly triggering pain.
- Loosening and migration of any or all of the components
- Alteration of neurological functions, occurrence of radiculopathies, dual lesions and/or pain, neurovascular deficit, including paralysis or other serious lesions. Effusion of cerebrospinal fluid.
- Digestive, urological and/or genital disorders, including stertion, impotence.
- Neurological injuries:

 - Hemorrhage and/or hematoma
 - Deep venous thrombosis, thrombophlebitis and/or pulmonary embolism
 - Adjacent segment disease
 - Disk inflammation, arachnoiditis and/or other types of inflammation.
 - Complications at bone graft donor site.
 - Inability to resume usual day-to-day activities.
 - Lesion of tissue issues.
 - Death

Note: Some potential adverse events may require an additional surgical revision procedure.

Storage and packaging:

The packaging of each of the components must be intact on delivery. A "tamper-proof" safety label testifies to the integrity of the packaging. If the word OPENED is visible, then the packaging must be considered as damaged, and it is not to be used.

The operating room personnel must thoroughly verify, before use, that all devices are complete and that no component presents any signs of damage. Damaged packaging and/or products must not be used but must be returned to SpineVision®.

The sterility use-by date must be complied with. Make sure that the implants are not damaged when opening the packaging. Damaged packages or products must not be used. If in doubt, please send them back to SpineVision®.

PRECAUTIONS:

Pedicle screw spinal systems must be implanted by an experienced spinal surgeon specifically trained in their use, as the technique required for their implantation is demanding, and exposes the patient to a serious risk of lesion. SURGICAL IMPLANTS MUST NEVER BE REUSED: An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage.

Preoperative:

- Implants are to be handled and stored with the utmost care. They must not be scratched or damaged. The implants and instruments must be protected while in storage, particularly against corrosive environments.
- The surgeon must be familiar with all the components before using the equipment, and must personally make sure that all the implants and instruments required are available before the start of the operation.
- The type of construction to be performed in each case must be determined before the start of the operation. An appropriate range of implants of various sizes, including smaller and larger than those initially provided for, must be available at the time of the operation.
- All components must be sterile before they can be used. Extra sterile components must be available, should they be needed unexpectedly.

Intraoperative:

- The surgeon must strictly follow the instructions for use of spinal instrumentation.
- The surgeon must be extremely careful to avoid any lesion to the spinal cord or to the nerve roots. This warning is of particular importance when inserting the screws. Any nervous lesion may induce a neurological dysfunction.
- A rupture, slippage or improper use of the instruments or implant components may cause injury to the patient or to members of the surgical team.
- The rigid parts of the Flex+2 rods must not be repeatedly or excessively bent beyond what is strictly necessary.
- The slightest stress on the Urethane Polycarbonate part of Flex+2 rods is prohibited, whether during bending of the rod or while the instrument is in place.
- In particular, make sure no pre-stressing is applied to the Urethane Polycarbonate part of the Flex+2 rods between the screws.
- The screw tightening components must not rest on the Urethane Polycarbonate part, as this might damage it.
- If the rods are to be kept to a certain length, the section must be a flat non-abrasive surface perpendicular to the axis. Cut the rods under the surgical field.
- Do not use screws with unsuitable dimensions (length, diameter), as this might damage the nerve roots, or cause hemorrhages and/or tearing.

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- Before closing, all locking screws must be tightened as per the instructions.

Postoperative:

- The postoperative recommendations given by the surgeon to the patient, and compliance with these recommendations by the patient, are extremely important.
- The patient must be informed of the limitations in the use of the device. He must be aware that complications such as deformation, coming apart or rupture of the device may occur if the spinal column is subjected to premature and/or excessive loads.
- The patient must be aware that motion of the spinal column will be limited at the instrumented level, and will have to be trained in making up for this permanent physical limitation in bodily motions.
- A persistent absence of bone setting at the fusion level which cannot be immobilized extends repeated and excessive stresses on the implant. Due to the material fatigue phenomenon, these stresses may cause the device to come apart, be distorted or ruptured. It is important to mobilize the fusion area and to confirm this setting through an X-ray examination. In case of abnormal absence of setting or of coming apart, distortion and/or rupture, the device is to be removed immediately, before it causes a serious lesion.

All explanted devices must be treated so as to make it impossible for them to be re-used during another surgical operation.

Magnetic Resonance (MRI) Safety:

The Flex+2 implants have not been evaluated for safety and compatibility in the MR environment. The Flex+2 implants have not been tested for heating or migration in the MR environment.

Disposal:

Disposal of all infectious wastes (implants) must be in accordance with applicable national and local regulations and the medical waste protocol of the medical center, and not endanger the security or the health of patient, user, or another person until its complete disposal (including cleaning, decontamination and waste).

Product Complaints:

Any health care professional (a.e., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identity, reliability, durability, safety, effectiveness and/or performance should notify SpineVision®. Further, if any of the implanted spinal instrumentation components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately. If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot, and date of use, your name and address, the nature of the complaint and notification as whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

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Email: corp.quality@spinevision.com

	See package insert for labeling limitation		Ethylene Oxide sterilisation
	Do not use if package is damaged and consult instructions for use		Consult instructions for use or consult electronic instructions for use
	Medical device		Batch code
	Do not re-use		Manufacturer
	Keep dry		Catalogue number
	Keep away from sunlight		Federal law (US) restricts this device to sale by or on the order of a physician

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