ENGLISH (OUT OF USA)

SpineVision®
Hexanium T1LF: Transforaminal Lumbar Intervertebral Body Fusion Device

Intended Use

The Hexanium T1LF (Transforaminal Lumbar Interbody Fusion) system is intended as intervertebral body fusion device of the lumbar spine (L2-S1). It is intended to be implanted via a transforaminal approach to give a primary stability of the anterior column (with or without restoring disc height) until fusion occurs.

Hexanium T1LF intervertebral body fusion device must be combined with an appropriate posterior lumbar fixation device, e.g. the SpineVision® UNI-Thread®, PLUS®, ULIS® or LUMIS®.

Description

The Hexanium T1LF system consists of a cage of various widths and heights, which can be inserted between two vertebral bodies of L2-S1 segment to give support and correction during lumbar interbody fusion surgeries. The hollow geometry of the implants allows them to be packed with bone graft material. This device has the potential of a primary stability of the anterior column (with or without restoring disc height) until fusion occurs.

For the Hexanium T1LF system, one cage must be used per segment for fusion in order to stabilize the segment completely.

The choice of graft will always be the surgeon's responsibility.

All cages are manufactured from implant grade Titanium (Ti6Al4V compliant with ASTM F3001).

According to the clinical evaluation report, it can be concluded that the benefit-risk ratio is positive.

Indications

The Hexanium T1LF (Transforaminal Lumbar Interbody Fusion) system is intervertebral body fusion device indicated for use with bone graft in skeletally mature patients with Degenerative Disc Disease (DDD) at one or two continuous levels from L2-S1. DDD is defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies. These DDD patients may also have up to Grade I spondylolisthesis or retrolisthesis at the involved level(s). Patients should have received at least 6 weeks of non-operative treatment prior to treatment with Hexanium T1LF system. This device has to be filled with bone graft material. This device has to be implanted via transforaminal approach. The Hexanium T1LF system must be used in combination with supplemental internal spinal fixation which has been cleared for use in the lumbar spine.

Contraindications

This device is not intended for thoracic or cervical spine use.

Contraindications include:

- infection, local to the operative site;
- signs of local inflammation;
- fever or leukocytosis;
- morbid obesity;
- pregnancy;
- pediatric cases, or patient still having general skeletal growth;
- spondylolisthesis unable to be reduced to Grade I;
- suspected or documented allergy or intolerance to composite metal;
- any case where the implant components selected for use would be too large or too small to achieve a successful result;
- any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality;
- any patient in which implant utilization would interfere with anatomical structures or expected physiological performance;
- bone fusion at the level to be treated;
- any case not needing a bone graft or fusion;
- any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis;
- any patient condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors or congenital abnormalities, fracture local to the operating site, elevation of segmentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count;
- mental illness;
- any patient unwilling to cooperate with postoperative instructions.

The contraindications of these devices are similar to those of other spinal interbody systems. This spinal instrument is not designed, or intended, or sold for uses other than those indicated. These contraindications must be taken into account by the physician when making his decision.

Potential Adverse Events

In addition to the risks associated with spinal surgery without instrumentation, the following potential adverse events can occur:

- delayed union;
- non-union (or pseudarthrosis);
- bone fracture or stress shielding at, above, or below the level of surgery, microfracture, resorption, or damage of any spinal bone (including pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site;
- retropulsed graft;
- loosening, or migration of implant, that may result in paralysis, neurological disorders or blood vessel erosion due to the proximity of the device;
- loss of neurological function, appearance of radiculopathy, dural tears, and/or development of pain, neurovascular compromise including paralysis, temporary or permanent retrograde ejaculation in males, or other types of serious injury, cerebral spinal fluid leakage
- infection;
- hemorrhage of blood vessels and/or hematomas;
- discitis, arachnoiditis, and/or other types of inflammation;
- subsidence of implant into adjacent vertebral bodies resulting in loss of height and possible impingement of neural structures;
- loss of proper spinal curvature, correction, height, and/or reduction;
- pain, discomfort, or other abnormal sensations due to the presence of the device;
- herniated nucleus pulposus, disc disruption or degeneration above, or below the level of surgery;
- loss of or increase in spinal mobility or function;
- pressure on the surrounding tissues or organs;
- foreign body reaction due to the presence of implants, such as a mass, auto-immune disease and/or impaired healing;
- cessation of any potential growth of the operated portion of the spine;
- deep venous thrombosis, thrombophlebitis;
- development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.;
- bone graft donor site complication;
- scar formation possibly causing neurological compromise or compression around nerves and/or pain;
- change in mental status;
- incapacity to resume the activities of normal everyday life;
- Death

Note: Some potential adverse events may require an additional surgical revision procedure.

Storage and packaging

Packages for each of the components should be intact upon receipt. Care should be taken not to damage the implants when opening the packaging. For implant supplied in a sterile form, carefully check that the seal on the cardboard box is intact before you open the box. Damaged packages or products must not be used. If in doubt, please send them back to SpineVision®.

Inspection, cleaning and sterilization

Hexanium T1LF implants are delivered sterile by exposure to a minimum dose of 25kGy of gamma radiation. They must be never resterilized.

Hexanium T1LF implants must only be used with the Hexanium T1LF ancillary equipment. Hexanium T1LF ancillary equipment are supplied not sterile. They must be sterilized before use.

All packaging and labeling materials must be removed prior to sterilization. Decontamination and cleaning must be performed before information regarding cleaning and decontamination of SpineVision® products are available in the instruction I-1G09 "Instructions for cleaning, sterilization, inspection and maintenance of SpineVision medical devices".

- Cleaning and decontamination of the non-sterile instruments at the Healthcare facility:** It is recommended to use a neutral pH, enzymatic or alkaline (pH < 11) cleaning agents solution and deionized or distilled warm (room temperature) water of soaking, cleaning and rinsing.

Disassemble the necessary instruments. Articulated instruments must be opened. Inspect each device and repeat the washing in the presence of any residual debris, paying close attention to threads, cannulas, hinges and hard to reach areas.

Note: Certain cleaning solutions, such as those containing bleach or formaldehyde, may damage certain components and must not be used.

Verification

The devices must always be verified before use: any devices presenting signs of weakness or surface scratches must not be used. (Return to the manufacturer or distributor any implants showing deterioration)

- Sterilization**

Prior to use, Hexanium T1LF instruments should be steam sterilized according to the following parameters:

Method	Cycle	Temperature	Minimum exposure time	Minimum drying time
Steam	Prevacuum	134° (273°F)	18 minutes	30 minutes

Validations were performed with dynamic-air-removal Steam sterilization process.

Warnings and precautions

PREOPERATIVE

- The use of the Hexanium T1LF has to be done using dedicated adequate instrumentation by a trained orthopedic surgeon or neurosurgeon familiar with these recommendations as well as with the operating techniques relating to these implants.

- The surgeon must be aware of not only the medical and surgical aspects of the implant, but also of the material limitations of titanium implants. The patient must be instructed in the limitations of the implant, especially with regard to the weight bearing and other body stresses on the device before a solid fusion occurs. The patient must be instructed that non-compliance with postoperative instructions could lead to failure of the device. The patient must be instructed that device failure could lead to the need for additional surgery to remove failed components.

Only patients that meet the criteria described in the indications should be selected and patients corresponding to the aforementioned contra-indications should not be selected. The following factors may influence when selecting patients for internal stabilization devices:

- foreign body sensitivity;
- certain degenerative diseases;
- senility, mental illness, alcoholism, or drug abuse;
- obesity;
- patient's occupation or activity level;
- smoking

- Additional sterile components should be available in case of any unexpected need.

INTRAOPERATIVE (Please see surgical technique GT13-S1, latest applicable version at SpineVision®)

- Correct handling of the implants is extremely important. The implants must not be transported in a way not to incur any damage. The implants must not be scratched and must not collide during handling.

- Care must be taken when picking up the Hexanium T1LF implant, filling it with bone graft and inserting it between the vertebral bodies. Alterations may produce defects in surface finish and may become the local point for eventual breakage of the implant or for decreased efficacy in withstanding post-operative migration (e.g. when teeth of implant are damaged).

Correct selection of the implant size is extremely important. The chance of adequate stabilization is increased by the selection of the proper size of the implant regarding the anatomy and conditions of the patient. Proper selection of the implant can minimize, but not eliminate risks.

- The Hexanium T1LF components have been designed for use with this system. Do not use the Hexanium T1LF with any components not specifically recommended, or with any components from another manufacturer.

- The Hexanium T1LF cannot be inserted via an anterior approach.

- The Hexanium T1LF cannot be used in cervical and thoracic spine.

- Implant should be inserted into the peripheral cortical bone of the endplates to minimize the risk of post-operative subsidence.

- Surgical implants must never be reused.** An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage

POSTOPERATIVE

- The fusion mass matures and becomes able to share load with the implant. **Adequately instruct the patient in the appropriate postoperative care.** The patient must be instructed to reduce stress on the implants in order to avoid clinical problems that may accompany failed fixation. The patient's ability and willingness to follow instructions are one of the most important aspects of successful bone healing. The patient must be instructed in the limitations of the implant, and that physical activity and full weight bearing may cause premature failure of internal stabilization devices by migrating, loosening or fracture. The patient must be instructed that an implant does not have the strength of normal healthy bone, and that device failure may occur if excessive loading is placed on the implant. An active, debilitated, or demented patient who is unable to use orthotic weight bearing devices may be particularly at risk.
- No internal stabilization device can withstand activity levels equal to those withstood by normal healthy segment. No implant can be expected to withstand unsupported weight-bearing stress indefinitely.
- Implants may fail when subjected to the increased loading associated with delayed union or non-union. Internal stabilization devices are intended as load-sharing devices until normal healing occurs. If normal healing does not occur or is delayed, the device may break due to fatigue. Patients should be informed of the risks of implant failure.

Magnetic Resonance (MR) Safety: The Hexanium T1LF has not been evaluated for safety and compatibility in the MR environment. The Hexanium T1LF has not been tested for heating or migration in the MR environment.

Disposal:

Disposal of all infectious wastes (implants) must be in accordance with applicable national and local regulations and the medical waste protocol of the medical center, and not endanger the security or the health of patient, user, or another person until its complete disposal (including cleaning, decontamination and waste).


Product complaints

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identification, durability, reliability, safety, effectiveness and/or performance, should notify SpineVision®. Further, if any of the implanted Hexanium T1LF components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately.

If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

In Europe	In USA
SpineVision® 10 Rue de la Renaissance 92160 Antony France Tél.: +33 (0) 1 5333 2525 Fax: +33 (0) 1 5333 2539 Email: corp.quality@spinevision.com	SpineVision Inc. 155 Federal Street 6th Floor – Suite 604 Boston, MA 02110- USA Phone: +1 888 770 4118 Fax: +1 888 770 4119

IMPORTANT INFORMATION CONCERNING HEXANIUM T1LF INSTRUMENTATION

In Europe	In USA
SpineVision® 10 Rue de la Renaissance 92160 Antony France Tél.: +33 (0) 1 5333 2525 Fax: +33 (0) 1 5333 2539 Email: corp.quality@spinevision.com	SpineVision Inc. 155 Federal Street 6th Floor – Suite 604 Boston, MA 02110- USA Phone: +1 888 770 4118 Fax: +1 888 770 4119

HEXANIUM T1LF: Transforaminal Lumbar Intervertebral Body Fusion Device

Purpose
The Hexanium T1LF (Transforaminal Lumbar Interbody Fusion) system is intended as intervertebral body fusion device of the lumbar spine (L2-S1). It is intended to be implanted via a transforaminal approach to give a primary stability of the anterior column (with or without restoring disc height) until fusion occurs.

Hexanium T1LF intervertebral body fusion device must be combined with an appropriate posterior lumbar fixation device, e.g. the SpineVision® UNI-Thread®, PLUS®, ULIS® or LUMIS®.

Description

The Hexanium T1LF system consists of a cage of various widths and heights, which can be inserted between two vertebral bodies of L2-S1 segment to give support and correction during lumbar interbody fusion surgeries. The hollow geometry of the implants allows them to be packed with autogenous bone graft.

For the Hexanium T1LF system, one cage must be used per segment for fusion in order to stabilize the segment completely.

The Hexanium T1LF must be filled with autogenous bone graft or a bone graft substitute to enhance the bone fusion. The choice of graft will always be the surgeon's responsibility.

All cages are manufactured from implant grade Titanium (Ti6Al4V compliant with ASTM F136).

Indications

The Hexanium T1LF (Transforaminal Lumbar Interbody Fusion) system is intervertebral body fusion device indicated for use with autogenous bone graft in skeletally mature patients with Degenerative Disc Disease (DDD) at one or two continuous levels from L2-S1. DDD is defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies. These DDD patients may also have up to Grade I spondylolisthesis or retrolisthesis at the involved level(s). Patients should have received at least 6 months of non-operative treatment prior to treatment with Hexanium T1LF system. This device has to be filled with autogenous bone graft material. This device is implanted via transforaminal approach.

Note for US audience only: The Hexanium T1LF system must be used in combination with supplemental internal spinal fixation which has been cleared by the FDA use in the lumbar spine

Contraindications

This device is not intended for thoracic or cervical spine use.

Contraindications include, but are not limited to:

- infection, local to the operative site;
- signs of local inflammation;
- fever or leukocytosis;
- morbid obesity;
- pregnancy;
- pediatric cases, or patient still having general skeletal growth;
- spondylolisthesis unable to be reduced to Grade I;
- suspected or documented allergy or intolerance to composite materials;
- any case where the implant components selected for use would be too large or too small to achieve a successful result;
- any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality;
- any patient in which implant utilization would interfere with anatomical structures or expected physiological performance;
- prior fusion at the level to be treated;
- any case not needing a bone graft or fusion;
- any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis;
- any other condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors or congenital abnormalities, fracture local to the operating site, elevation of segmentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count;
- mental illness;
- any patient unwilling to cooperate with postoperative instructions.

The contraindications of these devices are similar to those of other spinal interbody systems. This spinal instrument is not designed, or intended, or sold for uses other than those indicated.

Potential Adverse Events

In addition to the risks associated with spinal surgery without instrumentation, the following potential adverse events can occur (non-exhaustive list):

- delayed union;
- non-union (or pseudarthrosis);
- bone fracture or stress shielding at, above, or below the level of surgery, microfracture, resorption, or damage of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site;
- retropulsed graft;
- loosening, or migration of implant, that may result in paralysis, neurological disorders or blood vessel erosion due to the proximity of the device;
- loss of neurological function, appearance of radiculopathy, dural tears, and/or development of pain, neurovascular compromise including paralysis, temporary or permanent retrograde ejaculation in males, or other types of serious injury, cerebral spinal fluid leakage
- infection;
- hemorrhage of blood vessels and/or hematomas;
- discitis, arachnoiditis, and/or other types of inflammation;

- subsidence of implant into adjacent vertebral bodies resulting in loss of height and possible impingement of neural structures;
- loss of proper spinal curvature, correction, height, and/or reduction;
- pain, discomfort, or other abnormal sensations due to the presence of the device;
- herniated nucleus pulposus, disc disruption or degeneration above, or below the level of surgery;
- loss of or increase in spinal mobility or function;
- pressure on the surrounding tissues or organs;
- foreign body reaction due to the presence of implants, such as a mass, auto-immune disease and/or impaired healing;
- cessation of any potential growth of the operated portion of the spine;
- deep venous thrombosis, thrombophlebitis;
- development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.;
- bone graft donor site complication;
- scar formation possibly causing neurological compromise or compression around nerves and/or pain;
- change in mental status;
- incapacity to resume the activities of normal everyday life;
- death.

Note: Some potential adverse events may require an additional surgical revision procedure.

Storage and packaging

Packages for each of the components should be intact upon receipt. Care should be taken not to damage the implants when opening the packaging. For implant supplied in a sterile form, carefully check that the seal on the cardboard box is intact before you open the box. Damaged packages or products must not be used. If in doubt, please send them back to SpineVision®.

Cleaning and sterilization

Hexanium T1LF implants are delivered sterile by exposure to a minimum dose of 25kGy of gamma radiation. They must be never resterilized.

Hexanium T1LF implants must only be used with the Hexanium T1LF ancillary equipment. Hexanium T1LF ancillary equipment are supplied not sterile. They must be sterilized before use.

All packaging and labeling materials must be removed prior to sterilization. Decontamination and cleaning must be performed before sterilization. For details regarding the following sections, please refer to instructions for cleaning, sterilization, inspection and maintenance (Lubrication) of SpineVision® medical devices provided with the sets.

For cleaning and decontamination of the non-sterile instruments at the Healthcare facility:

It is recommended to use a neutral pH, enzymatic and alkaline (pH < 11) cleaning agents solution and deionized or distilled warm (room temperature) water of soaking, cleaning and rinsing.

Disassemble the necessary instruments. Articulated instruments must be opened. Inspect each device and repeat the washing in the presence of any residual debris, paying close attention to threads, cannulas, hinges and hard to reach areas.

Note: Certain cleaning solutions, such as those containing bleach or formaldehyde, may damage certain components and must not be used, except for aluminum instruments used on high-risk patients.

Note for European audiences only – New requirements related to high risk patients:

In the case of high risk patients, i.e. being suspected of contact with Non-conventional Transmissible Agents (e.g. anthrax), decontamination should be performed according to requirements of the World Health Organization (WHO). In this case, SpineVision® instruments can be decontaminated with sodium hydroxide, with the strict exception of instruments made out of aluminum, where bleach should be used.

Sterilization

Prior to use, Hexanium T1LF instruments should be steam sterilized according to the following parameters:

Method	Cycle	Temperature	exposure time	drying time
Steam	Prevacuum	134* (273°F) See note 1	18 minutes	30 minutes
Steam	Prevacuum	132* (270°F)	4 minutes	30 minutes

Validations were performed with dynamic-air-removal steam sterilization process.

Note 1: This sterilization cycle is not considered by the Food and Drug Administration to be a standard sterilization cycle. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have cleared by the Food and Drug Administration for the selected sterilization cycle specifications (time and temperature).

Lubrication specifications

SpineVision® recommends the lubrication instruments with moving or articulating parts. The instruments must be lubricated after each pre-sterilization / decontamination and before each steam sterilization.

Verification

The devices must always be verified before use: any devices presenting signs of weakness or surface scratches must not be used. (Return to the manufacturer or distributor any implants showing deterioration)

Warnings and precautions

PREOPERATIVE

- The use of the Hexanium T1LF has to be done using dedicated adequate instrumentation by a trained orthopedic surgeon or neurosurgeon familiar with these recommendations as well as with the operating techniques relating to these implants.
- The surgeon must be aware of not only the medical and surgical aspects of the implant, but also of the material limitations of titanium implants. The patient must be instructed in the limitations of the implant, especially with regard to the weight bearing and other body stresses on the device before a solid fusion occurs. The patient must be instructed that non-compliance with postoperative instructions could lead to failure of the device. The patient must be instructed that device failure could lead to the need for additional surgery to remove failed components.
- Only patients that meet the criteria described in the indications should be selected and patients corresponding to the aforementioned contra-indications should not be selected. The following factors may be of importance when selecting patients for internal stabilization devices:

- foreign body sensitivity;
- certain degenerative diseases;
- senility, mental illness, alcoholism, or drug abuse;
- obesity;
- patient's occupation or activity level;
- smoking

- Additional sterile components should be available in case of any unexpected need.

INTRAOPERATIVE

- Correct handling of the implants is extremely important. The implants must be transported in a way not to incur any damage. The implants must not be scratched and must not collide during handling.
- Care must be taken when picking up the Hexanium T1LF implant, filling it with autogenous bone graft or bone substitute and inserting it between the vertebral bodies. Alterations may produce defects in surface finish and may become the local point for eventual breakage of the implant or for decreased efficacy in withstanding post-operative migration (e.g. when teeth of implant are damaged).

- Correct selection of the implant size is extremely important. The chance of adequate stabilization is increased by the selection of the proper size of the implant regarding the anatomy and conditions of the patient. Proper selection of the implant can minimize, but not eliminate risks.

The Hexanium T1LF components have been designed for use with this system. Do not use the Hexanium T1LF with any components not specifically recommended, or with any components from another manufacturer.

- The Hexanium T1LF cannot be inserted via an anterior approach.

- The Hexanium T1LF cannot be used in cervical and thoracic spine.

- The implant should be placed on the peripheral cortical bone of the endplates to minimize the risk of post-operative subsidence.

- Surgical implants must never be reused.** An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage

POSTOPERATIVE

- Careful patient monitoring for the first two to three hours following the operation is very important while the fusion mass matures and becomes able to share load with the implant. **Adequately instruct the patient in the appropriate postoperative care.** The patient must be instructed to reduce stress on the implants in order to avoid clinical problems that may accompany failed fixation. The patient's ability and willingness to follow instructions are one of the most important aspects of successful bone healing. The patient must be instructed in the limitations of the implant, and that physical activity and full weight bearing may cause premature failure of internal stabilization devices by migrating, loosening or fracture. The patient must be instructed that an implant does not have the strength of normal healthy bone, and that device failure may occur if excessive loading is placed on the implant. An active, debilitated, or demented patient who is unable to use orthotic weight bearing devices may be particularly at risk.
- No internal stabilization device can withstand activity levels equal to those withstood by normal healthy segment. No implant can be expected to withstand unsupported weight-bearing stress indefinitely.
- Implants may fail when subjected to the increased loading associated with delayed union or non-union. Internal stabilization devices are intended as load-sharing devices until normal healing occurs. If normal healing does not occur or is delayed, the device may break due to fatigue. Patients should be informed of the risks of implant failure.

Magnetic Resonance (MR) Safety: The Hexanium T1LF has not been evaluated for safety and compatibility in the MR environment. The Hexanium T1LF has not been tested for heating or migration in the MR environment.


Product complaints

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identification, reliability, safety, effectiveness and/or performance, should notify SpineVision®. Further, if any of the implanted Hexanium T1LF components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately.

If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

In Europe	In USA
SpineVision® 10 Rue de la Renaissance 92160 Antony France Tél.: +33 (0) 1 5333 2525 Fax: +33 (0) 1 5333 2539 Email: corp.quality@spinevision.com	SpineVision Inc. 155 Federal Street 6th Floor – Suite 604 Boston, MA 02110- USA Phone: +1 888 770 4118 Fax: +1 888 770 4119

IMPORTANT INFORMATION CONCERNING HEXANIUM T1LF INSTRUMENTATION

In Europe	In USA
SpineVision® 10 Rue de la Renaissance 92160 Antony France Tél.: +33 (0) 1 5333 2525 Fax: +33 (0) 1 5333 2539 Email: corp.quality@spinevision.com	SpineVision Inc. 155 Federal Street 6th Floor – Suite 604 Boston, MA 0

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